

Der Spezialist für Naturheilmittel und Kosmetik. Seit **1903** im Dienste der Gesundheit.



Hyaluronic Acid Ampoules

The secret of our highly effective hyaluronic acid ampoules and their quickly visible "lifting effect" lies in the combination of short- and long-chain hyaluronic acid in conjunction with valuable moisturising and skin-firming active ingredients!

Long-chain (high-molecular) hyaluronic acid cannot penetrate the skin barrier and, together with the skin's keratin, forms a moisture-preserving film that simultaneously improves the skin's elasticity.

Short-chain (low-molecular) hyaluronic acid, on the other hand, has the ability to penetrate deeper into the skin and provide additional moisture in the long term. This plumping makes the skin appear smoother and firmer and, when used regularly, visibly reduces unsightly wrinkles and fine lines.

Aloe Vera gives the skin valuable moisture, protects and soothes irritated skin.

A **micro-algae** called **Nannochloropsis Oculata** is the basis for an active ingredient with visible skin-firming properties, which provides an immediate and long-lasting lifting effect. In combination with special polysaccharides, the formation of collagen-I is simultaneously stimulated; an indispensable component of the skin connective tissue, which protects and stabilises the skin.

A **special moisturising** factor based on plant sugar has been proven to improve skin hydration for many hours, contributing to the visibly smoothing effect that makes your skin look fresher, younger and more radiant.

The highly effective hyaluronic acid ampoules are suitable for daily use and can be dosed individually.

Ingredients: Aqua, Aloe Barbadensis Leaf Juice, Sorbitol, Saccharide Isomerate, Polysorbate 20, Sodium Hyaluronate, Aloe Barbadensis Leaf Extract, Nannochloropsis Oculata Extract, Pullulan, Benzoic Acid, Potassium Sorbate, Citric Acid, Sodium Citrate, Phenoxyethanol, Dehydroacetic Acid, Sodium Benzoate, Parfum.

Made in Germany

Contents 14 ampoules with 2 ml each for external use

Art.-Nr.: 1206

BP0520